


30" x 75"

76.2 cm x 190.5 cm


standard

12" x 12"
30.5 x 30.5 cm
window



premium

36" x 24"
91.4 x 70 cm
window




36" x 84"

91.4 cm x 213.4 cm


standard

12" x 12"
30.5 x 30.5 cm
window



premium

36" x 24"
91.4 x 70 cm
window




All mobile shielding dimensions are measured from the floor to the top of the paneling

Shielding available in 0.5 mm, 0.8 mm, or 1.6 mm lead equivalent (LE)

STRAIGHT BOOTH


standard

12" x 12"
30.5 x 30.5 cm
window



premium



48" x 24"
121.9 x 70 cm
window



ANGLE BOOTH


standard

12" x 12"
30.5 x 30.5 cm
window



premium

22.5" x 41.5"
57 x 105.4 cm
window



Shielding available in 0.5 mm, 0.8 mm, or 1.6mm lead equivalent (LE)