

VetCAT IQ™

High-Resolution CT Imaging for Bone & Soft Tissue

Main Specifications

Dimensions	32 W x 47 L X 60 / 67 H in 81W x 119 L x 152/171 H cm
Weight	472 lbs / 204 kg
Effective Bore Size	25 in / 63.5 cm diameter
Gantry Rotation	360 degrees

Power

Line Voltage	~120 V / ~230V
Line Frequency	50-60 Hz
Line Current	5 A standby / 2.5 A standby
Peak Current	20 A max / 10 A max

Acquisition Protocols

	Frames	Scan Time	3D Image Viewing
Bone - Large FOV	600	20 seconds	22 x 13 cm
Tissue - Large FOV	600	45 seconds	22 x 13 cm
Typical Reconstruction		10 - 76 seconds	
High-Resolution (default)		0.3 mm voxels in 3D	
Region of Interest (ROI)		≥ 70 microns	



Roll, Scan, and Diagnose in Seconds!

iM3 USA & Canada & Europe

US: info@im3usa.com | +1 800 664 6348



EU: sales@im3dental.com | support@im3dental.com

Manufactured by 

iM3® | THE GLOBAL NAME IN
VETERINARY DENTISTRY

With **VetCAT IQ**, you can see disease you might otherwise miss.

Scan image shown:
Tumor invading the orbital floor



iM3 USA & Canada & Europe

US: info@im3usa.com | +1 800 664 6348

EU: sales@im3dental.com | support@im3dental.com

Manufactured by **XORAN**

iM3 THE GLOBAL NAME IN
VETERINARY DENTISTRY