



the imaging solution for your clinic

## A Full Body High-Resolution CT Designed for All Veterinarians

### Main Specifications

Dimensions	45 W x 29 L X 65 H in 115 W x 74 L x 166 H cm
Weight	495 lbs / 224 kg
Effective Bore Size	31.5 in / 80 cm diameter
Gantry Rotation	360 degrees

### Power

Line Voltage	~120 / 240 V
Line Frequency	50 - 60 Hz
Line Current	5A / 2 A Standby
Peak Current	20 A / 10 A Max

### Acquisition Protocols

	Frames	Scan Time	FOV
Bone - large FOV	600	20 seconds	22 x 13 cm
Tissue - large FOV	600	45 seconds	22 x 13 cm

Typical reconstruction time	10 - 76 seconds
High-resolution (default)	0.3 mm voxels in 3D
Region of Interest (ROI)	≥ 70 microns (bone only)



Roll, Scan, and Diagnose in Seconds!

im3 USA & Canada & Europe

US: [info@im3usa.com](mailto:info@im3usa.com) | +1 800 664 6348


EU: [sales@im3dental.com](mailto:sales@im3dental.com) | [support@im3dental.com](mailto:support@im3dental.com)

Manufactured by XORAN


**im3** | THE GLOBAL NAME IN  
VETERINARY DENTISTRY



## Sample Scan Image



## Patient Positioning



## vTRON with sample patient bed



iM3 USA & Canada & Europe  
US: [info@im3usa.com](mailto:info@im3usa.com) | +1 800 664 6348  
EU: [sales@im3dental.com](mailto:sales@im3dental.com) | [support@im3dental.com](mailto:support@im3dental.com)

Manufactured by XORAN

**iM3**® | THE GLOBAL NAME IN  
VETERINARY DENTISTRY